

# 6-14-19 Initial practice for growth protocol

Thursday, June 6, 2019 10:27 AM

Start with 3 Pseudomonas isolates  
1 biological replicates  
2 technical replicates

## Day 1:

- ☐ 1. Streak isolates onto BHI agar plate
- ☐ 2. Incubate at 32°C overnight (16-18 h)

## Day 2:

- ☐ 1. Inoculate 5 mL BHI broth with a single colony
- ☐ 2. Vortex for 5 seconds
- ☐ 3. Incubate tube at 32°C overnight (16-18h)

## Day 3:

- ☐ 1. Vortex culture tube for 5 seconds
- ☐ 2. Perform serial dilutions in PBS (transferring 100 uL into 900 uL PBS) to the  $10^{-5}$  dilution
- ☐ 3. Spiral plate 50 uL of sample in duplicate onto BHI agar plates and incubate at 32°C overnight (16-18h)
- ☐ 4. Store inoculated BHI tubes at 4°C overnight (these will be used to inoculate milk on day 4).

## Day 4:

- ☐ 1. Enumerate cfu on BHI plates using Q count from your overnight culture and record in excel spreadsheet.
  - Past data showed that bacterial numbers ranged between  $7.0 \times 10^7$  and  $1.3 \times 10^8$  cfu/mL
- ☐ 2. Adjust with PBS to achieve an inoculum concentration of **100 cfu/mL**
- ☐ 3. Prepare inoculums: resuspend the desired culture in media (SMB milk). Transfer 1 mL of resuspended culture to 9 mL of uninoculated media (SMB milk), for a total of 3 tubes per media tested. Vortex for 5 seconds to resuspend.
- ☐ 4. Remove 500 uL for plating and, using the spiral plater, spiral plate both an undiluted culture (from step 2) and the  $10^{-1}$  dilution (from step 4) onto BHI plates in duplicate.
- ☐ 5. Incubate inoculated milk tubes at 6°C, and the spiral-plated agar plates at 32°C overnight (16-18 h). Be sure to include 1 uninoculated control tube which will be used to account for potential contaminants. (8AM)
- ☐ 6. Remove tubes from the 6°C incubator and vortex for 5 seconds. (8PM)
- ☐ 7. Aseptically transfer 500 uL of cultures into 1.5 mL Eppendorf tubes. Return tubes to 6°C incubator.
- ☐ 8. Perform serial dilutions to  $10^{-5}$ .
- ☐ 9. Plate undiluted,  $10^{-2}$ , and  $10^{-4}$  dilutions
- ☐ 10. Incubate inoculated plates at 32°C for 24 hrs.

## Day 5-6

- ☐ 11. Enumerate starting inoculum using the Q-count, to ensure the inoculum was **100 cfu/mL**. Record results in excel spreadsheet.
- ☐ 1. If the inoculum was between **90 – 200 cfu/mL** the tubes will continue to be incubated, if not, repeat days 2-5.
- ☐ 2. Remove tubes from the 6°C incubator and vortex for 5 seconds. (8AM)
- ☐ 3. Aseptically transfer 500 uL of cultures into 1.5 mL Eppendorf tubes. Return tubes to 6°C incubator.
- ☐ 4. Perform serial dilutions to  $10^{-5}$ .
- ☐ 5. Plate undiluted,  $10^{-2}$ , and  $10^{-4}$  dilutions
- ☐ 6. Incubate inoculated plates at 32°C for 24 hrs.
- ☐ 7. Remove tubes from the 6°C incubator and vortex for 5 seconds. (8PM)
- ☐ 8. Aseptically transfer 500 uL of cultures into 1.5 mL Eppendorf tubes. Return tubes to 6°C incubator.
- ☐ 9. Perform serial dilutions to  $10^{-5}$ .
- ☐ 10. Plate undiluted,  $10^{-2}$ , and  $10^{-4}$  dilutions
- ☐ 11. Incubate inoculated plates at 32°C for 24 hrs.

## Day 7

- ☐ 1. Enumerate inoculum using the Q-count. Record results in excel spreadsheet.
- ☐ 2. Remove tubes from the 6°C incubator and vortex for 5 seconds. (8AM)
- ☐ 3. Aseptically transfer 500 uL of cultures into 1.5 mL Eppendorf tubes. Return tubes to 6°C incubator.
- ☐ 4. Perform serial dilutions to  $10^{-5}$ .
- ☐ 5. Plate undiluted,  $10^{-2}$ , and  $10^{-4}$  dilutions
- ☐ 6. Incubate inoculated plates at 32°C for 24 hrs.

#### Day 8

- ☐ 7. Repeat day 7 until Nmax is reached
  - ☐ • When the count doesn't increase in two platings in a row
- ☐ 1. Calculate the growth rate by subtracting the initial-day concentration from the final-day concentration
- ☐ 2. Build growth curve